



Analysis steps

Beta bursts


1. Preprocess data
 - a. Reject bad electrodes
 - b. Bandpass filter (0.5-100Hz)
 - c. Notch filter (60 Hz)
 - d. Resample to 250Hz
 - e. Re-reference to an average reference
 - f. Remove bad segments through visual inspection
 - g. ICA for eye blinks and muscle artifacts
2. Import preprocessed raw file in database
 - a. Motor data – epoch time (-1.5s,10.5s)




b. Rest data – Split in blocks of 3.6s



3. Reject first 5 trials from rest and motor epochs
4. Concatenate time of raw files and extract channels of interest






5. Export the concatenated file to **.csv with header** following the proper naming convention:

```
# filename = participant_task_session_block.csv
# e.g. AD0109_motor_20Hz_post15.csv
# e.g. AD0109_rest_sham_base.csv

participant (e.g.AD0109), task (rest vs motor), session (20Hz vs 70Hz vs sham), block (base, post15, post45)
```



Brainstorm 24-Jan-2023

File Colormaps Update Plugins Help

George_Kenya_preprocessed

No data loaded.

Export matrix...

Save in: rest_3-6.sec

Recent It...

Desktop

Documents

This PC

Network

File name: AD0109_rest_20Hz_post15.csv

Files of type: ASCII: Comma-separated with header (*.csv)

Save

Cancel

Files to process: Data [0]

Process1 Process2

AutoSave: Off

AD0109_motor_70Hz_base... Saved to this PC

File Home Insert Draw Page Layout Formulas Data Review View Automate Help

Clipboard Font Alignment Number Styles Cells Editing Sensitivity Add-ins Analyze Data

POSSIBLE DATA LOSS Some features might be lost if you save this workbook in the comma-delimited (.csv) format. To preserve these features, save it in an Excel file format.

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S	T	U	V	W
1	Time	-1.5	-1.496	-1.492	-1.488	-1.484	-1.48	-1.476	-1.472	-1.468	-1.464	-1.46	-1.456	-1.452	-1.448	-1.444	-1.44	-1.436	-1.432	-1.428	-1.424	-1.42	-1.4
2	FC5	2.22E-06	-1.57E-06	-1.32E-06	-1.23E-06	-2.74E-06	-5.96E-07	1.21E-06	2.09E-06	3.11E-06	-7.99E-08	3.70E-07	2.92E-06	1.47E-06	2.84E-06	1.24E-06	-1.30E-06	6.94E-07	-2.21E-06	-2.35E-06	1.36E-07	-4.19E-06	-3.59E-
3	FC1	1.09E-07	-2.15E-06	-2.06E-06	-1.42E-06	-2.82E-06	-1.48E-06	1.31E-06	2.34E-06	1.97E-06	-7.85E-07	4.79E-07	3.48E-06	2.77E-06	3.39E-06	1.52E-06	-2.17E-07	1.16E-06	-2.14E-06	-2.87E-06	-1.03E-06	-3.77E-06	-2.66E-
4	C3	-6.57E-07	-3.39E-06	-3.09E-06	-4.46E-06	-3.64E-06	-4.82E-07	4.83E-07	8.72E-07	1.11E-06	-3.86E-07	8.76E-08	1.19E-06	1.12E-06	2.24E-06	4.90E-07	-1.27E-06	-2.29E-07	-2.57E-06	-3.43E-06	-2.05E-06	-4.14E-06	-3.67E-
5	CP5	1.37E-07	-3.23E-06	-2.92E-06	-2.83E-06	-4.02E-06	-1.78E-06	-1.08E-07	7.36E-07	1.63E-06	-1.41E-06	-5.46E-07	2.05E-06	5.91E-07	1.66E-06	2.11E-07	-1.67E-06	7.45E-07	-1.06E-06	-2.24E-06	-1.71E-06	-4.92E-06	-4.25E-
6	CP1	1.26E-08	-3.68E-06	-4.20E-06	-4.54E-06	-4.15E-06	-4.53E-07	6.07E-07	8.20E-07	8.02E-07	-2.26E-06	-1.98E-07	2.27E-06	8.14E-07	1.87E-06	-7.41E-07	-2.43E-06	7.12E-09	-3.27E-06	-4.02E-06	-2.07E-06	-5.21E-06	-4.36E-
7	FC3	-3.03E-06	-5.16E-06	-6.65E-06	-5.94E-06	-3.33E-06	-1.93E-07	3.01E-06	7.37E-06	7.89E-06	6.12E-06	6.28E-06	9.38E-06	7.83E-06	5.94E-06	1.88E-06	1.23E-07	4.47E-07	-2.25E-06	-2.53E-06	-1.87E-06	-5.04E-06	-3.14E-
8	C5	-2.53E-07	-5.13E-06	-3.76E-06	-1.77E-06	-1.17E-06	6.32E-07	1.50E-06	2.35E-06	5.73E-07	-2.28E-06	9.92E-07	2.58E-06	1.01E-06	1.96E-06	-5.97E-07	-1.16E-06	9.18E-07	-2.39E-06	-1.18E-06	4.01E-07	-4.51E-06	-3.18E-
9	C1	-4.64E-07	-4.27E-07	-3.29E-06	-7.89E-06	-6.72E-06	-3.83E-06	-1.24E-07	5.80E-06	6.17E-06	2.58E-06	2.68E-06	3.49E-06	2.85E-06	2.94E-06	1.95E-06	-4.25E-07	-1.35E-06	-1.93E-06	-2.15E-06	-1.28E-06	-2.33E-06	-2.20E-
10	CP3	3.56E-06	8.06E-07	9.95E-07	2.22E-06	1.11E-06	3.09E-06	5.91E-06	6.85E-06	6.20E-06	2.45E-06	1.26E-06	3.10E-06	2.41E-06	6.95E-07	8.57E-07	3.47E-06	5.31E-06	4.78E-06	4.58E-06	3.67E-06	1.37E-06	3.64E-

6. Run `extract_bursts_filtering_entire_trial.py` to extract bursts
 - a. This code will extract the burst characteristics (rate, amplitude, duration) for each of the above file and will be saved as `bursts_FILENAME.csv`
 - b. Before running, make sure to change the parameters as needed, in particular:
 - i. `import_folder` (where the files from step 5. were saved)
 - ii. `export_folder` (where the beta burst files will be saved)
 - iii. `mov_epoch_intervals` (depends on the motor task and how epochs were separated in Brainstorm)
 1. [-1.1s, -0.1s] Pre-movement
 2. [0.5s, 3.5s] Movement
 3. [5s, 8s] Post-movement
 - iv. `mov_epoch_name_intervals`
 - v. `rest_intervals` (depends on the resting state epochs were separated in Brainstorm)
 - vi. `rest_name_intervals`
 - c. The `bursts_FILENAME.csv` has each line representing one burst within the file. The first column is the burst identifier: `channel_interval_epoch_index`
 - i. Channel = channel from which the burst was extracted
 - ii. Interval = Movement or Resting state interval as defined in `mov_epoch_name_intervals` and `rest_name_intervals`
 - iii. Epoch = In which motor or rest epoch the burst is identified
 - iv. Index = The index number of the burst in the current channel, interval and epoch

AutoSave OFF | bursts_AD0109_motor_20Hz_post15.csv | Search

File Home Insert Draw Page Layout Formulas Data Review View Automate Help


Clipboard | Font: Aptos Narrow, 11 | Alignment: | Number: General

POSSIBLE DATA LOSS: Some features might be lost if you save this workbook in the comma-delimited (.csv) format. To preserve these features:

	A	B	C	D	E	F	G	H	I	J	K
1		rate	ampl	duration							
2	FC5_p-pre_1_0	0	0	0							
3	FC5_pre_1_0	0	0	0							
4	FC5_prep_1_0	0	0	0							
5	FC5_mov_1_0	0.666667	1.58E-06	148							
6	FC5_mov_1_1	0.666667	4.21E-06	116							
7	FC5_rel_1_0	0.666667	3.82E-07	112							
8	FC5_post_1_0	0.333333	8.35E-06	288							
9	FC5_p-post_1_0	0.798722	6.71E-06	192							
10	FC5_p-post_1_1	0.798722	2.80E-06	112							
11	FC5_p-pre_2_0	0	0	0							
12	FC5_pre_2_0	0	0	0							
13	FC5_prep_2_0	1.666667	7.40E-06	192							
14	FC5_mov_2_0	0.333333	7.58E-07	140							
15	FC5_rel_2_0	0.666667	4.88E-07	116							
16	FC5_post_2_0	0.333333	4.59E-06	116							
17	FC5_p-post_2_0	0.399361	1.15E-05	328							
18	FC5_p-pre_3_0	0	0	0							
19	FC5_pre_3_0	0	0	0							
20	FC5_prep_3_0	0	0	0							
21	FC5_mov_3_0	0.333333	1.70E-06	120							
22	FC5_rel_3_0	0	0	0							
23	FC5_post_3_0	0	0	0							
24	FC5_p-post_3_0	1.198083	3.05E-06	112							
25	FC5_p-post_3_1	1.198083	4.66E-06	204							
26	FC5 p-post 3 2	1.198083	3.25E-06	112							

7. Create the bad_channels.csv file

- Run `bad_channels.py` and set the “path” variable as the export folder used above.
- Open the created file in Excel and manually enter on the second column all the bad channels.
 - Bad channels can be obtained in Brainstorm by right clicking on an epoch from each file > View all the bad channels
 - If no bad channels are within a file, write “N”



	A	B	C	D	E	F	G
1	BM1212_rest_tDCS_post15.csv	F8, TP9, TP10, PO9, PO10					
2	BM1212_motor_tDCS_post15.csv	TP9, TP10, PO9, PO10					
3	BM1212_motor_tDCS_post45.csv	TP9, TP10, PO9, PO10, AF8					
4	BM1212_rest_tDCS_post45.csv	F8, TP9, TP10, PO9, PO10					
5	BM1212_motor_tDCS_base.csv	F8, AF8					
6	BM1212_rest_tDCS_base.csv	N					
7	BM1212_rest_sham_post15.csv	Fz, F4, F8, FC5					
8	BM1212_motor_sham_post15.csv	Fz, F4, F8, FC5, FC2					
9	BM1212_rest_sham_post45.csv	Fz, F8					
10	BM1212_motor_sham_post45.csv	Fz, F4, FC5					
11	BM1212_rest_sham_base.csv	Fz, F4, F8, FC5					
12	BM1212_motor_sham_base.csv	Fz, F4, FC5					
13	BM1212_rest_tACS_post15.csv	Fz, F4, FC5, FC1, C3, CP1					
14	BM1212_motor_tACS_post15.csv	Fz, F4, FC5, TP10, PO10					
15	BM1212_rest_tACS_post45.csv	Fz, FC5					
16	BM1212_motor_tACS_post45.csv	F4, FC5					
17	BM1212_rest_tACS_base.csv	F4, F8, FC5					
18	BM1212_motor_tACS_base.csv	Fz, F4, FC5, TP10, PO10					
19	JD0213_motor_sham_post15.csv	FC1, TP9, TP10, O2, FC4					
20	JD0213_rest_sham_post15.csv	FC1, O1					
21	JD0213_motor_sham_post45.csv	FC1, Cz, TP10, FC4, CP4					
22	JD0213_rest_sham_post45.csv	FC1, O1, FC4					
23	JD0213_motor_sham_base.csv	N					
24	JD0213_rest_sham_base.csv	FT9					
25	JD0213_motor_tDCS_post15.csv	Fz, F4, F8, FC5, TP10, PO10					

8. Run `graph_bursts.py` to plot the bursts

- Note: `graph_bursts.py` needs to be changed based on what you need to plot
- Requires:
 - A `bad_channels.csv` file
 - `functions.py` file to be in the same folder
- Change parameters at the top of the file:
 - channels
 - motor_intervals_of_interest
 - rest_intervals_of_interest
- Comment/Uncomment the `line_error()` function calls to plot the data you are interested in